

Tammy Rellar

Americas Sales Manager Tekmar, Leeman, CETAC

Global Product Marketing Manager CETAC

Automation & Liquids





- Aerospace and Defense
- Factory Automation
- Air and Water Quality Environmental Monitoring
- Oceanographic Research

- Energy
- Medical Imaging
- Pharmaceutical Research
- Electronics Design



• Driving Technology: Improving our Night Vision

 Tracking Air Quality Onboard an Aircraft in Real Time

Improving Water Safety and Early Virus Detection

Measuring Ocean Conditions in Storms

Detecting Cancer and Treating It

Temperature Screening by Thermal Imaging

Tracking Ocean Garbage from Space

Understanding Asteroids

Teledyne Technologies

Instrumentation Digital Imaging Aerospace and Defense Electronics Engineered Systems

Instrumentation


Our Instrumentation segment provides measurement, monitoring and control instruments for marine, environmental, scientific and industrial applications. We also provide power and communications connectivity devices for distributed instrumentation systems and sensor networks deployed in mission critical, harsh environments.



- 25 Businesses
- Accounts for about 45% of overall Teledyne Technologies sales
- Presented as 4 market segments: Instruments, Water Quality, Marine and Oil & Gas
- Tekmar participates in the Instruments and Water Quality segments
- Tekmar is part of the EEMI Group Environmental and Electronic Measurement Instruments
 - Test & Measurement, ISCO, API, TAMS and TACS

Teledyne Environmental & Electronic Measurement Instruments (EEMI)

Vicki Benne- Chief Executive Officer



T&M: Teledyne Test and Measurement

ISCO: Teledyne ISCO

TACS: Teledyne Advanced Chemistry Systems

TAMS: Teledyne Advanced Monitoring Solutions

TAPI: Teledyne Advanced Pollution Instruments



Laboratory Instrumentation for Organic and Inorganic Analysis of Food, Water and Soil Quality:

- Sample Preparation and Concentration of Trace VOCs
- Total Organic Carbon (T)C)
- Measurements for Various Water Analyses

Instrumentation for Elemental Analysis

- Inductively Coupled Plasma Spectrometers
- DC Arc Spectrometers
- Mercury Analyzers

Instrumentation and Sensors for Measurement & Control of Gas Flows & Vacuum Pressures

- Thermal, Pirani, Piezo, & Ionization
 Vacuum Sensors
- Thermal Flow Sensors & Proportional Control valves

Accessories for Sample Introduction and Elemental Analysis

- Autosamplers
- Sample Preparation
- Nebulizers
- Laser Ablation Systems

OEM Automation Services

Laser-based Sample Introduction
Systems and Analytical Instruments

- Excimer Laser Ablation Systems
- Femtosecond Laser Ablation System
- Diode Laser Systems
- CO2 Laser Systems







Current Laboratory Challenges



- Pipette Tips not available
- Staffing Shortages
- Social Distancing

All have negative impact on sample preparation





What's so bad about sample preparation?







For the Analyst

- Repetitive Motion
 - Decapping
 - Transferring
 - Pipetting
 - Risk of Carpal Tunnel
- Seems to take all day
- Low job satisfaction
 - Leading to high turnover
- Boring!



IS THIS WHY I BECAME A CHEMIST???









Who made the standards?

For the Manager

- Nature of repetitive work leads to mistakes
 - Which then requires more sample preparation
- Time Consuming
- Even the best analysts:
 - Get Sick
 - Take Leave
- Paying someone to stand in one place all day...
- High Staff Turnover





Opportunity for Error....and why not use your analysts for more important things?



Introducing the SimPrep



- Based on CETAC automation
- Recognized excellence of Hamilton micro station
- Use same sample racks ICP autosamplers and SimPrep
- Simple easy to use software
- Small mobile footprint
- Full sample homogenization
- Aqueous sample handling including
 - Soil digests
 - Mining
 - Drinking Water
 - Waste Waters
 - List goes on....






Video of Operation







More than Just Dilution




- Dispensing
- Dosing
- Serial Dilutions
- Split samples
- Spike samples
- Prepare Standards
- And More




Precision SimPrep vs Hand Prep





6 Point Calibration Curve





- Hand Prepped via serial dilution technique
- Diluted from 100ppb stock
- Both Hand prep and SimPrep produce good calibration data
- On the hand preparation the curve was contaminated for Zinc

One Autosampler Two Flow Paths

- Switch the Autosampler between ICP and SimPrep software on the same PC
- (SimPrep) Automated calibration standards preparation
- (SimPrep) Automated Sample Dilution
- (Autosampler) Automated Sample Introduction
- Requires USB Hub





System Configuration Options







SimPrep - 280





Having trouble sourcing pipette tips? Just tired of paying for them?

- The SimPrep precision syringes remove the need for costly and possibly difficult to obtain pipette tips.
- Resulting in lower consumable costs and less plastic waste
- Automating the liquid transfer step removes an opportunity for error





Summary & Benefits of the SimPrep

- Automate tedious sample preparation steps
- Reduce consumable usage
- Free up staff for more critical tasks
- One autosampler for sample preparation and introduction
- Plug and play installation
- Intuitive software
- Unify and improve sample and standard preparation





Thank You!

For more information, please visit

https://www.teledynecetac.com/products/automation/simprep

